

CAB4429

Product Information

Product SKU: CAB4429	Gene ID: 10514	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 149kDa	Conjugate: Unconjugated
Calculated MW: 149kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a nucleolar transcriptional regulator that was first identified by its ability to bind specifically to the Myb proto-oncogene protein. The encoded protein is thought to play a role in many cellular processes including response to nucleolar stress, tumor suppression and synthesis of ribosomal DNA. Alternate splicing results in multiple transcript variants.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: P160; PAP2; Pol5; MYBBP1A

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1104-1328 of human MYBBP1A (NP_055335.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.