

## CAB4459

---

### Product Information

<b>Product SKU:</b> CAB4459	<b>Gene ID:</b> 10890	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Mouse

---

### Additional Information

<b>Observed MW:</b> 23kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 23kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** RAB10 belongs to the RAS (see HRAS; MIM 190020) superfamily of small GTPases. RAB proteins localize to exocytic and endocytic compartments and regulate intracellular vesicle trafficking (Bao et al., 1998 [PubMed 9918381]).

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** RAB10

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 98-198 of human RAB10 (NP\_057215.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.