

CAB4472

Product Information

Product SKU: CAB4472	Gene ID: 10991	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 60kDa	Conjugate: Unconjugated
Calculated MW: 56kDa	Isotype: IgG

Immunogen Information

Background: Enables L-amino acid transmembrane transporter activity. Involved in amino acid transport. Located in plasma membrane.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: G17; SN1; NAT1; SNAT3; DEE102; SLC38A3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-65 of human SLC38A3 (NP_006832.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.