

CAB4505

Product Information

Product SKU: CAB4505	Gene ID: 11333	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 26kDa/25kDa	Conjugate: Unconjugated
Calculated MW: 21kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is a phosphoprotein that may upregulate the PDGFA-stimulated growth of fibroblasts and also downregulate the mitogenicity of PDGFB. The encoded protein in rodents has been shown to bind PDGFA with a low affinity.

Recommended Dilution: WB,1:100 - 1:500 IF/ICC,1:50 - 1:200

Synonyms: PAP; PAP1; HASPP28; PDAP1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PDAP1 (NP_055706.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.