

CAB4513

Product Information

Product SKU: CAB4513	Gene ID: 22829	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 100kDa	Conjugate: Unconjugated
Calculated MW: 92kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a type I membrane protein that belongs to the family of neuroligins, which are cell adhesion molecules present at the postsynaptic side of the synapse, and may be essential for the formation of functional synapses. Alternatively spliced transcript variants have been found for this gene.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:100

Synonyms: HNL4Y; NLGN4Y

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 122-256 of human NLGN4Y (NP_001157710.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.