

CAB4528

Product Information

Product SKU:	CAB4528	Gene ID:	23122	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	141kDa	Conjugate:	Unconjugated
Calculated MW:	141kDa	Isotype:	IgG

Immunogen Information

Background: Enables cytoskeletal protein binding activity; dystroglycan binding activity; and protein tyrosine kinase binding activity. Involved in several processes, including microtubule cytoskeleton organization; positive regulation of extracellular matrix organization; and regulation of supramolecular fiber organization. Located in several cellular components, including basal cortex; cortical microtubule plus-end; and ruffle membrane. Colocalizes with focal adhesion; kinetochore; and microtubule cytoskeleton.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: CLASP2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 220-400 of human CLASP2 (NP_001193973.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.