

## CAB4540

---

### Product Information

<b>Product SKU:</b> CAB4540	<b>Gene ID:</b> 23283	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human,Mouse

---

### Additional Information

<b>Observed MW:</b> 70kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 64kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** Enables RNA binding activity. Predicted to be involved in pre-mRNA cleavage required for polyadenylation. Located in nucleoplasm.

**Recommended Dilution:** WB,1:1000 - 1:2000 IF/ICC,1:50 - 1:200

**Synonyms:** CstF-64T; CSTF2T

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 420-560 of human CSTF2T (NP\_056050.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.