

CAB4580

Product Information

Product SKU: CAB4580	Gene ID: 23761	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human, Mouse

Additional Information

Observed MW: 50kDa	Conjugate: Unconjugated
Calculated MW: 47kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene catalyzes the conversion of phosphatidylserine to phosphatidylethanolamine in the inner mitochondrial membrane. The encoded protein is active in phospholipid metabolism and interorganelle trafficking of phosphatidylserine.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: PSD; LIBF; PSDC; PSSC; DJ858B16; dJ858B16.2; PISD

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-375 of human PISD (NP_055153.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.