

CAB4581

Product Information

Product SKU: CAB4581	Gene ID: 24139	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human, Mouse

Additional Information

Observed MW: 71kDa	Conjugate: Unconjugated
Calculated MW: 71kDa	Isotype: IgG

Immunogen Information

Background: Enables microtubule binding activity. Involved in negative regulation of microtubule polymerization and regulation of microtubule nucleation. Colocalizes with mitotic spindle.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: ELP70; EMAP2; EMAP-2; EML2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 150-415 of human EML2 (NP_036287.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.