

CAB4600

Product Information

Product SKU:	CAB4600	Gene ID:	26002	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	70kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable copper ion binding activity and dopamine beta-monoxygenase activity. Predicted to be involved in dopamine catabolic process; norepinephrine biosynthetic process; and octopamine biosynthetic process. Part of endoplasmic reticulum membrane.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: MOX; PRO5780; dJ248E1.1; MOXD1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 390-590 of human MOXD1 (NP_056344.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.