

CAB4620

Product Information

Product SKU:	CAB4620	Gene ID:	26258	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	20kDa	Conjugate:	Unconjugated
Calculated MW:	20kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene may play a role in intracellular vesicle trafficking. It interacts with Syntaxin 13 which mediates intracellular membrane fusion. Mutations in this gene cause symptoms associated with Hermansky-Pudlak syndrome-9. Alternative splicing results in multiple transcript variants. A pseudogene related to this gene is located on the X chromosome.

Recommended Dilution: WB, 1:200 - 1:2000

Synonyms: PA; HPS9; PLDN; BLOS6; PALLID; BLOC1S6

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human BLOC1S6 (NP_036520.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.