

CAB4644

Product Information

Product SKU: CAB4644	Gene ID: 27125	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human, Mouse

Additional Information

Observed MW: 150kDa	Conjugate: Unconjugated
Calculated MW: 127kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene belongs to the AF4 family of transcription factors involved in leukemia. It is a component of the positive transcription elongation factor b (P-TEFb) complex. A chromosomal translocation involving this gene and MLL gene on chromosome 11 is found in infant acute lymphoblastic leukemia with ins(5;11)(q31;q31q23).

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: MCEF; CHOPS; AF5Q31; AFF4

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human AFF4 (NP_055238.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.