

CAB4655

Product Information

Product SKU: CAB4655	Gene ID: 3656	Size: 20uL, 100uL
Clone No: ARC1078	Host Species: Rabbit	Reactivity: Human,Rat

Additional Information

Observed MW: 70kDa	Conjugate: Unconjugated
Calculated MW: 69kDa	Isotype: IgG

Immunogen Information

Background: IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: IRAK-2; IRAK2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 50-250 of human IRAK2 (O43187).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.