

CAB4662

Product Information

Product SKU:	CAB4662	Gene ID:	28981	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	79kDa	Conjugate:	Unconjugated
Calculated MW:	80kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene, together with IFT74, forms a tubulin-binding module of intraflagellar transport complex B. This module is involved in transport of tubulin within the cilium, and the encoded protein is required for ciliogenesis. Mutations in this gene are a cause of short-rib polydactyly syndromes.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: DV1; CDV1; CDV-1; CDV1R; CDV-1R; SRTD19; IFT81

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 550-676 of human IFT81 (NP_054774.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.