

CAB4674

Product Information

Product SKU: CAB4674	Gene ID: 29796	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 7kDa	Conjugate: Unconjugated
Calculated MW: 7kDa	Isotype: IgG

Immunogen Information

Background: UCRC is a subunit of mitochondrial complex III (ubiquinol-cytochrome c reductase; EC 1.10.2.2), which forms the middle segment of the respiratory chain of the inner mitochondrial membrane (Schagger et al., 1995 [PubMed 8592474]).

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: QCR9; UCRC; HSPC051; HSPC119; HSPC151; UCCR7.2; UQCR10

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-63 of human UQCR10 (NP_037519.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.