

## CAB4694

---

### Product Information

<b>Product SKU:</b>	CAB4694	<b>Gene ID:</b>	50808	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	26kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	26kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The protein encoded by this gene is a GTP:ATP phosphotransferase that is found in the mitochondrial matrix. Several transcript variants encoding a few different isoforms have been found for this gene.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** AK6; FIX; AK3L1; AKL3L; AKL3L1; AK3

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-227 of human AK3 (NP\_057366.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.