

CAB4734

Product Information

Product SKU: CAB4734	Gene ID: 51388	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human, Mouse

Additional Information

Observed MW: 20kDa	Conjugate: Unconjugated
Calculated MW: 20kDa	Isotype: IgG

Immunogen Information

Background: Enables RNA binding activity. Predicted to be involved in ribosomal large subunit biogenesis. Located in cytosol; nucleolus; and nucleoplasm.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: KD93; CGI-37; HSPC031; NIP7

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human NIP7 (NP_057185.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.