

CAB4785

Product Information

Product SKU: CAB4785	Gene ID: 54507	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 116kDa	Conjugate: Unconjugated
Calculated MW: 117kDa	Isotype: IgG

Immunogen Information

Background: This gene is a member of ADAMTS (a disintegrin and metalloproteinase with thrombospondin motifs)-like gene family and encodes a protein with seven thrombospondin type 1 repeats. The thrombospondin type 1 repeat domain is found in many proteins with diverse biological functions including cellular adhesion, angiogenesis, and patterning of the developing nervous system. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: TSRC1; ECTOL2; ADAMTSL-4; ADAMTSL4

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 795-1074 of human ADAMTSL4 (NP_061905.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.