

CAB4842

Product Information

Product SKU: CAB4842	Gene ID: 55720	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human, Mouse

Additional Information

Observed MW: 110kDa	Conjugate: Unconjugated
Calculated MW: 92kDa	Isotype: IgG

Immunogen Information

Background: Enables RNA binding activity. Predicted to be involved in endonucleolytic cleavage of tricistronic rRNA transcript (SSU-rRNA, 5.8S rRNA, LSU-rRNA) and maturation of SSU-rRNA from tricistronic rRNA transcript (SSU-rRNA, 5.8S rRNA, LSU-rRNA). Predicted to be located in nucleolus. Predicted to be part of preribosome, small subunit precursor.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: TSR1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 455-804 of human TSR1 (NP_060598.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.