

CAB4853

Product Information

Product SKU: CAB4853	Gene ID: 55922	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 78kDa	Conjugate: Unconjugated
Calculated MW: 78kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a transcriptional repressor that interacts with specific negative regulatory elements to mediate transcriptional repression of certain nuclear factor kappa B responsive genes. The protein localizes predominantly to the nucleolus with a small fraction found in the nucleoplasm and cytoplasm. Alternate splicing results in multiple transcript variants.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: NRF; ITBA4; NKRF

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 391-690 of human NKRF (NP_060014.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.