

## CAB4864

---

### Product Information

<b>Product SKU:</b>	CAB4864	<b>Gene ID:</b>	56648	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	17kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	17kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to enable translation elongation factor activity. Predicted to be involved in positive regulation of translational elongation. Located in intracellular membrane-bounded organelle.

**Recommended Dilution:** WB, 1:200 - 1:2000

**Synonyms:** EIF-5A2; eIF5AII; EIF5A2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-153 of human EIF5A2 (NP\_065123.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.