

## CAB4897

---

### Product Information

<b>Product SKU:</b>	CAB4897	<b>Gene ID:</b>	57705	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	354kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	354kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to be involved in autophagy. Predicted to act upstream of or within with a positive effect on CD8-positive, alpha-beta T cell activation. Predicted to act upstream of or within antigen processing and presentation and cellular response to virus. Predicted to be located in early endosome and endoplasmic reticulum.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** C10orf64; WDFY4

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human WDFY4 (NP\_065996.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.