

CAB4902

Product Information

Product SKU: CAB4902	Gene ID: 57863	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 43kDa	Conjugate: Unconjugated
Calculated MW: 43kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is a calcium-independent cell-cell adhesion protein that can form homodimers or heterodimers with other nectin proteins. The encoded protein has both homophilic and heterophilic cell-cell adhesion activity. This gene is reported to be a tumor suppressor gene.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: BlgR; NECL1; TSLL1; CMT2FF; IGSF4B; Nectl-1; synCAM3; CADM3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-330 of human CADM3 (NP_001120645.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.