

CAB4905

Product Information

Product SKU: CAB4905	Gene ID: 58498	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Mouse

Additional Information

Observed MW: 19kDa	Conjugate: Unconjugated
Calculated MW: 19kDa	Isotype: IgG

Immunogen Information

Background: Predicted to enable calcium ion binding activity. Located in A band.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: MYL2A; MYLC2A; MYL7

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of human MYL7 (NP_067046.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.