

CAB4912

Product Information

| | | |
|-----------------------------|------------------------|--------------------------|
| Product SKU: CAB4912 | Gene ID: - | Size: 20uL, 100uL |
| Clone No: - | Host Species: - | Reactivity: - |

Additional Information

| | |
|-------------------------|---------------------|
| Observed MW: - | Conjugate: - |
| Calculated MW: - | Isotype: - |

Immunogen Information

Background: -

Recommended Dilution: -

Synonyms: -

Purification Method: -

Immunogen: -

Storage: -