

FBRS Rabbit Polyclonal Antibody



CAB4927

Product Information

Size:

50uL, 100uL

Observed MW:

53kDa

Calculated MW:

48kDa/51kDa

Applications:

WB IHC IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000 IHC 1:50
- 1:200 IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Fibrosin is a lymphokine secreted by activated lymphocytes that induces fibroblast proliferation (Prakash and Robbins, 1998 [PubMed 9809749]).[supplied by OMIM, Mar 2008]

Immunogen information

Gene ID:

64319

Uniprot

Q9HAH7

Synonyms:

FBRS; FBS; FBS1; fibrosin


Immunogen:

Recombinant protein of human FBRS

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of Jurkat cells, using FBRS antibody (CAB4927). Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.