

CAB4945

Product Information

Product SKU: CAB4945	Gene ID: 65003	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 21kDa	Conjugate: Unconjugated
Calculated MW: 21kDa	Isotype: IgG

Immunogen Information

Background: This nuclear gene encodes a 39S subunit component of the mitochondrial ribosome. Alternative splicing results in multiple transcript variants. Pseudogenes for this gene are found on chromosomes 5 and 12.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: L11MT; CGI-113; MRP-L11; MRPL11

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-192 of human MRPL11 (NP_057134.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.