

CAB4990

Product Information

| | | | | | |
|---------------------|---------|----------------------|--------|--------------------|-----------------|
| Product SKU: | CAB4990 | Gene ID: | 79982 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|-----------------------|-------|-------------------|--------------|
| Observed MW: | 42kDa | Conjugate: | Unconjugated |
| Calculated MW: | 43kDa | Isotype: | IgG |

Immunogen Information

Background: Enables Hsp70 protein binding activity. Involved in cellular protein-containing complex assembly and chaperone cofactor-dependent protein refolding. Located in endoplasmic reticulum.

Recommended Dilution: WB,1:1000 - 1:4000

Synonyms: EGNR9427; PRO34683; DNAJB14

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human DNAJB14 (NP_001026893.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.