

CAB5000

Product Information

Product SKU: CAB5000	Gene ID: 80279	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 69kDa	Conjugate: Unconjugated
Calculated MW: 57kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a protein that has been reported to function in signaling pathways governing transcriptional regulation and cell cycle progression. It may play a role in tumorigenesis and metastasis. A pseudogene of this gene is located on the long arm of chromosome 20. Alternative splicing results in multiple transcript variants that encode different isoforms.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:100

Synonyms: C53; IC53; LZAP; HSF-27; MST016; PP1553; OK/SW-cl.114; CDK5RAP3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CDK5RAP3 (NP_788276.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.