

CAB5017

Product Information

Product SKU: CAB5017	Gene ID: 81617	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human, Mouse

Additional Information

Observed MW: 37-39kDa	Conjugate: Unconjugated
Calculated MW: 39kDa	Isotype: IgG

Immunogen Information

Background: Predicted to enable protein serine/threonine kinase activator activity. Predicted to be involved in intracellular signal transduction. Predicted to be located in cytosol.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: MO2L; MO25-BETA; bA103J18.3; CAB39L

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-337 of human CAB39L (NP_112187.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.