

CAB5057

Product Information

Product SKU:	CAB5057	Gene ID:	8411	Size:	20uL, 100uL
Clone No:	ARC1250	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	170kDa	Conjugate:	Unconjugated
Calculated MW:	162kDa	Isotype:	IgG

Immunogen Information

Background: Enables 1-phosphatidylinositol binding activity; GTP-dependent protein binding activity; and protein homodimerization activity. Involved in endocytosis; vesicle fusion; and viral RNA genome replication. Located in cytosol and early endosome. Is extrinsic component of plasma membrane.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: MST105; ZFYVE2; MSTP105; EEA1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human EEA1 (Q15075).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.