

CAB5068

Product Information

Product SKU:	CAB5068	Gene ID:	6596	Size:	20uL, 100uL
Clone No:	ARC1982	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	114kDa	Conjugate:	Unconjugated
Calculated MW:	114kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the SWI/SNF family. Members of this family have helicase and ATPase activities and are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein contains a RING finger DNA binding motif. Two transcript variants encoding the same protein have been found for this gene. However, use of an alternative translation start site produces an isoform that is truncated at the N-terminus compared to the full-length protein.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: ZBU1; HLTF1; RNF80; HIP116; SNF2L3; HIP116A; SMARCA3; HLTF

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 900-1009 of human HLTF (Q14527).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.