

CAB5080

Product Information

Product SKU: CAB5080	Gene ID: 5879/5880/5881	Size: 20uL, 100uL
Clone No: ARC1165	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 20kDa	Conjugate: Unconjugated
Calculated MW: 21kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: MIG5; MRD48; Rac-1; TC-25; p21-Rac1; Rac1/2/3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-189 of human RAC1(P63000).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.