

FUT11 Rabbit Polyclonal Antibody



CAB5190

Product Information

Size:

50uL

Observed MW:

Calculated MW:

53kDa/55kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

170384

Uniprot

Q495W5

Synonyms:

FUT11


Immunogen:

Recombinant protein of human FUT11

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using FUT11 antibody (CAB5190).