

CAB5375

Product Information

Product SKU: CAB5375	Gene ID: 3399	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: Refer to figures	Conjugate: Unconjugated
Calculated MW: 13kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is a helix-loop-helix (HLH) protein that can form heterodimers with other HLH proteins. However, the encoded protein lacks a basic DNA-binding domain and therefore inhibits the DNA binding of any HLH protein with which it interacts.

Recommended Dilution: IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:100

Synonyms: HEIR-1; bHLHb25; ID3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human ID3 (NP_002158.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.