

## CAB5586

---

### Product Information

<b>Product SKU:</b>	CAB5586	<b>Gene ID:</b>	1749	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC1994	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	32kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	32kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a member of a homeobox transcription factor gene family similar to the Drosophila distal-less gene. The encoded protein may play a role in bone development and fracture healing. Mutation in this gene, which is located in a tail-to-tail configuration with another member of the family on the long arm of chromosome 7, may be associated with split-hand/split-foot malformation.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** SHFM1; SHFM1D; DLX5

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 150-250 of human DLX5 (P56178).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.