

## CAB5930

---

### Product Information

<b>Product SKU:</b>	CAB5930	<b>Gene ID:</b>	51286	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2092	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	22kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	15kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The protein encoded by this gene is a neuron-specific protein. The similar protein in pig enhances neuroblastoma cell differentiation in vitro and may be involved in neuronal differentiation in vivo. Multiple pseudogenes have been reported for this gene.

**Recommended Dilution:** WB,1:1000 - 1:5000 IF/ICC,1:50 - 1:200

**Synonyms:** BM88; CEND1

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 70-149 of human CEND1 (NP\_057648.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.