

CAB5955

Product Information

Product SKU:	CAB5955	Gene ID:	51330	Size:	20uL, 100uL
Clone No:	ARC2119	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	17kDa	Conjugate:	Unconjugated
Calculated MW:	14kDa	Isotype:	IgG

Immunogen Information

Background: Involved in positive regulation of extrinsic apoptotic signaling pathway and regulation of wound healing. Predicted to be located in cell surface and ruffle. Predicted to be active in plasma membrane.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: FN14; CD266; TWEAKR; TWEAKR/Fn14/CD266

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 50-129 of human TWEAKR/Fn14/CD266 (Q9NP84).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.