

CAB5996

Product Information

Product SKU:	CAB5996	Gene ID:	6810	Size:	20uL, 100uL
Clone No:	ARC2113	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	34kDa	Conjugate:	Unconjugated
Calculated MW:	34kDa	Isotype:	IgG

Immunogen Information

Background: Enables sphingomyelin phosphodiesterase activator activity. Involved in several processes, including cornified envelope assembly; positive regulation of immune effector process; and positive regulation of protein localization. Located in several cellular components, including basolateral plasma membrane; cytoplasmic vesicle; and lamellipodium. Part of SNARE complex. Is active in glutamatergic synapse and postsynapse.

Recommended Dilution: WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

Synonyms: STX4A; p35-2; Syntaxin 4

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Syntaxin 4 (NP_004595.2)).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.