

CAB6004

Product Information

| | | |
|-----------------------------|-----------------------------|---------------------------------|
| Product SKU: CAB6004 | Gene ID: 124540 | Size: 20uL, 100uL |
| Clone No: - | Host Species: Rabbit | Reactivity: Human, Mouse |

Additional Information

| | |
|-----------------------------|--------------------------------|
| Observed MW: 35kDa | Conjugate: Unconjugated |
| Calculated MW: 35kDa | Isotype: IgG |

Immunogen Information

Background: This gene encodes an RNA-binding protein that is a member of the Musashi protein family. The encoded protein is transcriptional regulator that targets genes involved in development and cell cycle regulation. Mutations in this gene are associated with poor prognosis in certain types of cancers. This gene has also been shown to be rearranged in certain cancer cells.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: MSI2H; Musashi-2 (MSI2)

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human Musashi-2 (Musashi-2 (MSI2)) (NP_620412.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.