

CAB6174

Product Information

Product SKU:	CAB6174	Gene ID:	4673	Size:	20uL, 100uL
Clone No:	ARC1888	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	52kDa	Conjugate:	Unconjugated
Calculated MW:	45kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the nucleosome assembly protein (NAP) family. This protein participates in DNA replication and may play a role in modulating chromatin formation and contribute to the regulation of cell proliferation. Alternative splicing results in multiple transcript variants encoding different isoforms; however, not all have been fully described.

Recommended Dilution: WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

Synonyms: NRP; NAP1; NAP1L; NAP1L1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NAP1L1 (P55209).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.