

CAB7399

Product Information

Product SKU: CAB7399	Gene ID: 121457	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Mouse,Rat

Additional Information

Observed MW: 39kDa	Conjugate: Unconjugated
Calculated MW: 39kDa	Isotype: IgG

Immunogen Information

Background: Predicted to be involved in response to X-ray. Located in endoplasmic reticulum.

Recommended Dilution: WB,1:1000 - 1:5000 IF/ICC,1:100 - 1:500

Synonyms: IKIP; IKBIP

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 115-350 of human IKBIP (NP_963906.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.