

CAB8250

Product Information

Product SKU: CAB8250	Gene ID: 64342	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 55kDa	Conjugate: Unconjugated
Calculated MW: 43kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene shares similarity with mouse Hs1bp3, an Hcls1/Hs1-interacting protein that may be involved in lymphocyte activation.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: ETM2; HS1-BP3; HS1BP3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human HS1BP3 (NP_071905.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.