

## CAB8317

---

### Product Information

<b>Product SKU:</b> CAB8317	<b>Gene ID:</b> 80817	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human, Mouse

---

### Additional Information

<b>Observed MW:</b> 45kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 44kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** Enables microtubule binding activity. Involved in centriole replication and centriole-centriole cohesion. Located in centriole; centrosome; and spindle pole.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** PS1TP3; KIAA1712; CEP44

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human CEP44 (NP\_001035247.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.