

CAB8655

Product Information

Product SKU:	CAB8655	Gene ID:	6663	Size:	20uL, 100uL
Clone No:	ARC1767	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	62kDa	Conjugate:	Unconjugated
Calculated MW:	50kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional activator after forming a protein complex with other proteins. This protein acts as a nucleocytoplasmic shuttle protein and is important for neural crest and peripheral nervous system development. Mutations in this gene are associated with Waardenburg-Shah and Waardenburg-Hirschsprung disease.

Recommended Dilution: WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

Synonyms: DOM; WS4; PCWH; WS2E; WS4C; SOX-10; SOX10

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 350-450 of human SOX10 (P56693).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.