

## CAB8662

---

### Product Information

<b>Product SKU:</b>	CAB8662	<b>Gene ID:</b>	6675	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC1769	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	59kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	59kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables identical protein binding activity. Predicted to be involved in UDP-N-acetylglucosamine biosynthetic process. Located in cytosol; nucleoplasm; and plasma membrane.

**Recommended Dilution:** WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

**Synonyms:** AGX; AGX1; AGX2; SPAG2; UAP1

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 400-500 of human UAP1 (Q16222).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.