

## CAB8694

---

### Product Information

<b>Product SKU:</b>	CAB8694	<b>Gene ID:</b>	4985	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC1786	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	45kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	40kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables G protein-coupled enkephalin receptor activity. Involved in several processes, including G protein-coupled opioid receptor signaling pathway; cellular response to hypoxia; and positive regulation of peptidyl-serine phosphorylation. Is intrinsic component of plasma membrane.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** DOP; DOR; DOR1; OPRD; OPRD1

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 1-100 of human OPRD1 (P41143).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.