

CAB8708

Product Information

Product SKU: CAB8708	Gene ID: 3493	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: 60kDa	Conjugate: Unconjugated
Calculated MW: 38kDa	Isotype: IgG

Immunogen Information

Background: Contributes to immunoglobulin receptor binding activity. Involved in antibacterial humoral response; glomerular filtration; and positive regulation of respiratory burst. Located in extracellular space. Part of monomeric IgA immunoglobulin complex and secretory dimeric IgA immunoglobulin complex.

Recommended Dilution: WB,1:500 - 1:1000 FC,1:20 - 1:50

Synonyms: IgA1; IGHA1

Purification Method: Affinity purification

Immunogen: A synthetic Peptide of human IGHA1

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.