

CAB9091

Product Information

Product SKU:	CAB9091	Gene ID:	126393	Size:	20uL, 100uL
Clone No:	ARC1787	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	17kDa	Conjugate:	Unconjugated
Calculated MW:	17kDa	Isotype:	IgG

Immunogen Information

Background: This locus encodes a heat shock protein. The encoded protein likely plays a role in smooth muscle relaxation.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: HEL55; Hsp20; PPP1R91; HSP20/HSPB6

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 70-150 of human HSP20/HSPB6 (O14558).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.