

CAB9156

Product Information

Product SKU:	CAB9156	Gene ID:	7529/7534	Size:	20uL, 100uL
Clone No:	ARC1455	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	28kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a protein belonging to the 14-3-3 family of proteins, members of which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals. The encoded protein has been shown to interact with RAF1 and CDC25 phosphatases, suggesting that it may play a role in linking mitogenic signaling and the cell cycle machinery. Two transcript variants, which encode the same protein, have been identified for this gene. [provided by RefSeq, Jul 2008]

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: GW128; HEL-S-1; HS1; KCIP-1; YWHAA; 14-3-3-beta/zeta

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human 14-3-3 beta/zeta (P31946).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.